



## ADSORBABLE ORGANIC HALIDES (AOX) IN WATER

### INTRODUCTION

#### AOX – WHAT ARE THEY?

Halogenated organic compounds (HOC) have increasingly been used within various industries over the years. Some man-made halogenated compounds are toxic, mutagenic, or carcinogenic and they may have harmful effects on human health and the environment.

Adsorbable Organic Halogen (AOX) Compounds (X = Cl, Br, I) is a sum parameter for describing the organic halogen compound load in water, sewage sludge and soils. The parameter covers a large group of substances from simple volatile substances such as trichloromethane (chloroform), or complex organic molecules such as dioxins/furans with a large variety of toxic properties. Most AOXs are chlorine-containing molecules, but bromo- and iodo-AOXs do also occur.

AOX compounds pose a potential concern because they resist breaking down in the environment therefore they have long half-life periods. Hence, AOX can be important in effluent quality monitoring from landfill leachates or industry in order to meet discharge limits.

#### METHOD INFORMATION

Test Parameter: Adsorbable Organic Halogen (AOX)

Test Method: In House Method based on USEPA 1650

Limit of Reporting: 0.1 mg/L

Sample size required: 500mL

Turnaround time: 5-7 working days from the receipt of samples.

### SOURCES AND FORMATION OF AOX

One of the main sources of AOX discharges world-wide is the pulp and paper industry. When using chlorine and chlorine chemicals to bleach fibres, some of the chlorine reacts with organic matter to form chlorinated organic substances which may then be discharged with the treated effluent and detected as AOXs.

Landfills can also be a source of necessary precursor organic compounds. Landfill leachates rich in organics which are discharged into waterways may flow into the water supply system. Leachates containing aromatic compounds will easily react with active chlorine and may generate halogenated compounds.

Smaller quantities of AOX are also formed during the routine chlorination of drinking water, swimming pools, cooling waters and process waters in laundries. In water treatment, the use of disinfectants such as chlorine, calcium hypochlorite and sodium hypochlorite results in the formation of free chlorine and combined chlorine residuals and disinfection by-products. Although many specific chlorine disinfection by-products have been identified, several of the total organic halogens have yet to be identified.

Other sources of AOX discharges are the effluents of industrial facilities which contain chlorine-containing substances from manufacturing or use. Raw hospital wastewater and effluent is also considered as an important source of AOX in municipal wastewater [3]. The chemical industry and other sectors such as metal surface treatment, textiles, waste treatment, printing and dry cleaning are also potential sources.

Follow us:



[www.alsglobal.com](http://www.alsglobal.com)  
RIGHT SOLUTIONS RIGHT PARTNER

HQ Shah Alam +603 7845 8257  
Other Sites:  
Johor Bahru    Pulau Pinang  
Kota Kinabalu    Kemaman

## CLASSES OF HOC's

There are a series of different forms of halogenated organic compounds which can be determined:

### AOX (Adsorbable Organically Bound Halogens)

Represents the organically bound halogens - chlorine, bromine and iodine (but not fluorine) contained in a sample which can be adsorbed on activated carbon.

### EOX (Extractable Organically Bound Halogens)

Represents the organically bound halogens - chlorine, bromine and iodine (but not fluorine) - which can be extracted with a non-polar solvent.

### POX (Purgeable Organically Bound Halogens)

Represents the organically bound halogens - chlorine, bromine and iodine (but not fluorine) - contained in a sample, which can be purged in the gas phase under defined conditions with an auxiliary gas.

AOX is by far the most common of the classes which are monitored and analysed however this is highly dependent upon specific industries.

## ANALYSIS OF ADSORBABLE ORGANICS HALOGENATED COMPOUNDS

The method of determining AOX has been standardized worldwide through a large number of standards. The ALS methodology for the analysis of AOX is based on the US EPA Method 1650.

AOX/TOX methodology uses carbon adsorption with a microcoulometric-titration detector to detect all organic halides containing chlorine, bromine & iodine that are adsorbed by granular activated carbon under specific conditions.

The result is a parameter used for water quality control purposes and represents the sum of organically bound chlorine, bromine and iodine (but not fluorine) that can be adsorbed on activated carbon under specified conditions and, if the sample is not filtered, includes that associated with suspended matter.



## IMPACTS ON HUMAN HEALTH AND ENVIRONMENT

There is no direct relationship between AOX value and toxicity. Excessive exposure to chemicals in this diverse grouping of chemicals may affect health with the possible effects depending on the particular chemical. Many of the chemicals detected under the parameter AOX are also individually listed. Some chlorinated chemicals detected by AOX are toxic to aquatic organisms at low concentrations. Many are persistent and can also bio-accumulate.

AOX compounds pose a potential concern because they resist breaking down in the environment. Some of these molecules are toxic at high concentrations. Because they can accumulate in the food chain, they pose a potential threat to aquatic organisms living in estuaries near bleached pulp effluents.

## GENERAL SAMPLING & PRESERVATION REQUIREMENTS

Water sample: 500mL Amber Glass Bottle with Teflon Cap, Zero Headspace. pH<2 with HNO<sub>3</sub>.

Storage: Cool to 4°C

Holding Time: 6 months

Please contact us to discuss your analysis requirement and we will assist you with the proposal.

WhatsApp: 017-555 2985 (Abdul Qaiyum)

Email: [marketing.my@alsglobal.com](mailto:marketing.my@alsglobal.com)